

DEPiLiGHT

808nm Diode Laser

Remove your concern about unwanted hair and expect immediate skin tightening

Permanent hair removal & Skin tightening

○ Why DEPILIGHT?

- Reduces burning risk
- Easy to use
- Speedy and easy operation
- Versatile indications

○ What are the benefits?

- 1 You can treat as fast as you can.
DEPiLiGHT has high coverage rate with a large spot size 9x9mm and high repetition rate up to 10Hz (for FDP mode), up to 3Hz (for DP mode).
- 2 High power continuous wave diode.
Energy delivery homogeneous along pulse segments gives reliable results. Blend mode pulses enable safer and more effective treatments on all skin types.
- 3 Sapphire tip contact cooling.
Gives patient comforts despite of delivery of higher fluences.
- 4 Continuous wave technology.
it's much more advanced technology than Quasi CW diode lasers enables both very short and very long pulses.
- 5 360W peak power is the same as the average power, which provides optimal therapeutic range and no side effects.
- 6 No consumables.
- 7 Versatile indications.
Not only for hair removal, but also it shows immediate results on skin tightening and acne treatment.



GTG Wellness

www.gtgwellness.co.kr

DEPiLiGHT

808nm Diode Laser

Permanent hair removal
& Skin tightening

DEPiLiGHT Specifications

| | |
|------------------------|--|
| Light Source | Diode(Continuous Wave) |
| Wavelength | 808nm |
| Fluence | Up tp 60J/cm2 |
| Pulse Duration | 10 - 430ms |
| Spot Size | 9 x 9 mm |
| Repetition Rate | DP mode - Up to 3Hz FDP mode - Up to 10Hz |
| Delivery | Sapphire tip contact |
| Cooling Tip | Size : 20 x 19 mm 0°C Chiller cooling |
| Weight | 26kg (cart excluded) |
| Volume | 450 x 500 x 275 |
| Electrical Requirement | 100~230 V (Free Voltage) |

Clinical Data



APPLICATION

- Permanent hair removal
- Acne Scar and treatment
- Skin tightening
- Wrinkle care



Handpiece



LCD Screen

